



**Data Sheet: CellGro® / CellGenix™**  
**Recombinant Human**  
**Interleukin 10 (IL-10)**  
**research grade - Order No.: 1414**

Source	A DNA sequence encoding the mature IL-10 protein was expressed in <i>E. coli</i>
Raw Material Statement	n.a.
Production	Produced in a GMP environment
Molecular Mass	The methionyl form of the <i>E. coli</i> expressed mature human IL-10 contains 167 amino acid residues (including a C-terminal His6-tag) and has a molecular mass of 19.6 kDa
Purity	≥ 95%, as determined by SDS-PAGE
Endotoxin Level	< 0.1 ng per 1 µg (< 1000 EU/mg), as determined by LAL gel clot test
Activity	Measured in a cell proliferation assay using an IL-10 sensitive cell line, MC/9.  ED50 is typically 0.05 – 0.2 ng/ml. Activity: ≥2.5x10 <sup>6</sup> IU/mg, calibrated with the 1 <sup>st</sup> International Standard NIBSC, # 93/722
Sterility	Sterility test of the vialled product (direct inoculation)
Formulation	Lyophilized from a 0.2 µm -filtered solution in PBS
Amount per Vial	50 µg or 10 µg, as determined by spectrophotometrical measurement at 280 nm using PBS as reference
Transport	Ambient temperature
Storage at	-20°C to -80°C. Avoid repeated freeze-thaw cycles
Shelf life	2 years from date of production
Intended Use	<b>For research use only</b> <b>Not intended for diagnostic or therapeutic use</b>

Manufacturer  
CellGenix Technologie Transfer GmbH  
Am Flughafen 16 – 79108 Freiburg – Germany  
Tel: +49 761 88889 – 0  
Fax: +49 761 88889 – 830  
www.cellgenix.com - info@cellgenix.com